

Better chemistry. Better business.

Hubbard-Hall Inc. • 563 South Leonard Street, Waterbury, CT 07608 • HubbardHall.com • 800-648-3412

Product Bulletin

Product Name: Metal Guard 830 Product Code: 2508030 Revision Date: June 22, 2020

Metal Guard® 830

Metal Guard 830 is a non-foaming, water-soluble, corrosion inhibitor that readily forms an effective rust preventative for iron and steel alloys that are packaged and stored indoors. Metal Guard 830 provides excellent short-term rust protection on cast iron and powdered metal substrates

Metal Guard 830 does not contain any inorganic salts; therefore, no powder residues are left on the surface of the metal. The thin film left behind also burns off cleanly in annealing or other high temperature operations.

Features & Benefits

Effective at low concentrations	Reduces operational cost
Low foam	Can be used in spray wash systems

Physical Data

Specific gravity	1.02
Solubility in water	Infinite
Appearance and odor	Clear amber liquid
pH (conc.)	8.19

Typical Applications

In-process short-term corrosion inhibitor

Operating Conditions

Metal Guard 830 can be applied by immersion, spray, or flow-coat. Drying can be accelerated by heating the Metal Guard 830 solution or by drying in a heated atmosphere.

Concentration	1% – 10% (volume)
Temperature	Ambient – 160° F (71°C)









Better chemistry. Better business.

Hubbard-Hall Inc. • 563 South Leonard Street, Waterbury, CT 07608 • HubbardHall.com • 800-648-3412

Product Bulletin

Product Name: Metal Guard 830 Product Code: 2508030 Revision Date: June 22, 2020

Corrosion Protection Properties

Metal Guard 830 may be classified as a non-severe inhibitor, that is, it should not be used on parts or fabrications that are stored outdoors under exposure to the elements and temperature extremes. Metal Guard 830 solutions should be used on parts for in-process protection stored indoors or boxed after the Metal Guard 830 solution has dried.

Titration Method

- 1. Pipette a 25 mL sample into a 250 mL Erlenmeyer flask.
- 2. Add 20 drops of 0.04% Methyl Red indicator.
- 3. Titrate with 0.1 N Hydrochloric Acid until color changes from yellow to red.
- 4. Record mL used.

Calculation

Concentration = mL 1.0 N HCl x 1.01

Test Kit Method

- 1. Using syringe, measure a 5 mL sample of the bath and dispense into the mixing bottle.
- 2. Add about 20 mL water and 5 to 10 drops Methyl Orange indicator.
- 3. Add 0.72 N Hydrochloric Acid dropwise until the solution turns pink.
- 4. Record the number of drops used.

Calculation

Concentration = # Drops 0.72 N HCl x 0.71

Waste Disposal

Discharge to a disposal system. In order to be completely informed on the latest regulations for your area, please contact the local authorities.









Better chemistry. Better business.

Hubbard-Hall Inc. • 563 South Leonard Street, Waterbury, CT 07608 • HubbardHall.com • 800-648-3412

Product Bulletin

Product Name: Metal Guard 830 Product Code: 2508030 Revision Date: June 22, 2020

WARRANTY: THE QUALITY OF THIS PRODUCT IS GUARANTEED ON SHIPMENT FROM OUR PLANT. IF THE USE RECOMMENDATIONS ARE FOLLOWED, DESIRED RESULTS WILL BE OBTAINED. SINCE THE USE OF OUR PRODUCTS IS BEYOND OUR CONTROL, NO GUARANTEE EXPRESSED OR IMPLIED IS MADE AS TO THE EFFECTS OF SUCH USE, OR THE RESULTS TO BE OBTAINED.

Our people. Your problem solvers.

For more information on this process please call us at 1-800-648-3412

or techservice@hubbardhall.com











